CE PSP Single Use Nitrile Glove
 Model: PSP 50-225, PSP 50-230
 Cat III Chemische en biologische risico's

President Safety B.V.,
 Postbus 100, 3220 AC Hellevoetsluis / NL
 www.pspasafety.com

A. GEBRUIK

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Déclaration de conformité UE : pspasafety.com/nl/declaration-of-conformity

Version nr.: 03-2021 ▲

CE PSP Single Use Nitrile Glove
 Code: PSP 50-225, PSP 50-230
 Category III - Chemical and biological risks

President Safety B.V.,
 Postbus 100, 3220 AC Hellevoetsluis / NL
 www.pspasafety.com

A. USE

The glove types named here meet with the requirements of regulation (EU) 2016/425, EN ISO 374-1:2016 and EN 420:2003+A1:2009. As the gloves are adjusted to special purposes, the lengths may deviate from the requirements of EN 420:2003+A1:2009. 100 Pieces (by weight) packed in dispenser.

RESULTS CHEMICAL TESTS:

EN ISO 374-1 TYPE B		EN ISO 374-5				
RESULTS CHEMICAL TESTS:						
n-heptane (J): Level 3 40% Sodium Hydroxide (K): Level 6 37% Formaldehyde (T): Level 6						
Level	1	2	3	4	5	6
Breakthrough time (min)	> 10	> 30	> 60	> 120	> 240	> 480
EN 374-4:2013 - Resistant to degradation against chemicals:						
n-heptane (J): 38,8 % (+/- 4,8 %) 40% Sodium Hydroxide (K): -6,8 % (+/- 4,4 %) 37% formaldehyde (T): 23,7 % (+/- 9,9 %)						
PPE is subject to the conformity assessment procedure Module D under surveillance of the notified body: SGS Fimko Oy, P.O. Box 30 (Sarkiniementie 3), 00211 Helsinki, Finland. Notified Body no: 0598						
EU type examination carried out by: CENTEXBEL-GENT, Technologiepark 70, 9052 Zwijnaarde, BELGIUM. NB. nr: 0493						

SIZING ACCORDING TO EN 420:2003+A1:2009

Handsize Index		Size				
	Label	07	08	09	10	11
		S	M	L	XL	XXL

B. WARNING / RISK ASSESSMENT:
 a) Gloves to protect against mechanical action whose effects are superficial, substances and mixtures which are hazardous to health, and harmful biological agents. Important: The gloves are recommended for use in situations where only low chemical protection is identified as needed. b) While selecting an equipment, user should perform risk analysis based on the intended use and determine the suitability based on product's test standards and protection levels obtained. c) Information provide does not reflect the actual duration of protection in the workplace due to other factors influencing the performance, such as temperature, abrasion, and degradation, and the differentiation between mixtures and pure chemicals d) Information regarding protection refers to the working surface, i.e. the palm of the glove, which has been submitted to testing. e) Gloves should be thoroughly inspected for damages (specially for nicks and holes) before use. If any damage is found avoid usage. f) The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm - where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture. g) It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation h) When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves i) The maximum wear time depends on the activity being carried out and the person. j) EN374-4:2013 Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemical. k) The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen. l) This product contains Nitrile Butadiene Rubber and compounding chemicals which may cause allergic reaction in some individuals. m) The glove offers no protection against perforation with sharp objects, e.g. injection needles.

STORAGE: To be stored in cool dry place away from sunlight.Storage procedures are the main factor in determining glove shelf life. Gloves should be kept in their packaging protected from sunlight, artificial light, humidity and stored at temperatures between 5 °C - 30 °C. Real time stability testing for expiration date on this products after accelerated ageing process are declared as 5 years.

Further information and declaration:
 President Safety B.V.,
 P.O. Box 100, 3220 AC Hellevoetsluis / The Netherlands
 www.pspasafety.com

CE REGULATION (EU) 2016/425
 0598

Déclaration de conformité : pspasafety.com/nl/declaration-of-conformity

Version nr.: 03-2021 ▲